DCOOK

IRRIGADOR BUCAL PORTÁTIL

IRRIGADOR BUCAL PORTÁTIL

PORTABLE ORAL IRRIGATOR

IRRIGATEUR BUCAL PORTABLE

MANUAL DE INSTRUCCIONES

MANUAL DE INSTRUÇÕES • INSTRUCTION MANUAL
NOTICE D'UTILISATION • MANUALE DI ISTRUZIONI

IRRIGADOR BUCAL PORTÁTIL

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PORTABLE ORAL IRRIGATOR

DESCRIPTION

1. Operating mode selector, smooth, normal and strong
2. On/off button.
3. Mode of operation indicator light
4. Water tank
5. USB charging cable
6. 5 heads
 - a. 2x Standard Nozzle
 - b. 1x Tongue Cleaner Nozzle (Tongue Scraper)
 - c. 1x Periodontal Tip
 - d. 1x Orthodontic Tip

SAFETY WARNINGS

• Please read these instructions carefully before using the product and keep them for future reference.

• This device is not intended for use by children under the age of 8, persons with reduced physical, sensory, or mental capabilities, or inexperience, unless they are supervised or have received appropriate instructions.

• Do not use the irrigator if you have wounds in the mouth, bleeding gums without a diagnosed cause or recent surgery in the oral area, unless medically indicated.

• Consult your dentist before using this product if you have dental implants, orthodontics, periodontal disease, or have had recent dental work.

• Do not aim the water directly at the eyes, ears, nasal cavities, or at other people or animals.

• Do not insert any parts of the device into body orifices other than the mouth.

- Use only drinking water or solutions recommended by the manufacturer. Do not add rinses with alcohol, essential oils, bleach, or other unauthorized products.
- Do not use the device while charging. Charge only with the adapter or cable supplied by the manufacturer.
- Do not immerse the device in water or expose it to direct jets. Although it is waterproof, it is not submersible.
- Do not use if the appliance, nozzle or charging cable is damaged. If in doubt, contact an authorized service center.
- Do not disassemble or attempt to repair the product yourself. Repairs should be carried out only by qualified personnel.
- Keep the device out of the reach of children when not in use. Small nozzles can pose a choking hazard.
- Use only compatible nozzles recommended by the manufacturer. The use of non-genuine accessories may affect the operation and safety of the device.

- Discontinue use and consult a professional if you experience pain, excessive bleeding, or unusual discomfort.
- This device is intended for personal and household use only. Do not use for medical or professional purposes without proper approval.
- Keep the plug and appliance dry, do not handle it with wet hands and always unplug it before cleaning it.
- Use, charge and store the appliance in dry, protected environments at temperatures between 15°C and 35°C. Do not expose it to the sun or cover it during charging.
- Use only the USB charging cable provided in this device.

BEFORE FIRST USE

- Remove all packaging and protections.
- Wash the water tank and nozzle with warm water before first use.
- Fully charge the device using the supplied USB cable. The full charge time is usually 3 hours. Check the LED indicator for charging status.

Note: Do not use the irrigator while it is charging.

HOW TO USE

1. Open the water tank lid.
2. Fill the reservoir with drinking water. You can use warm water for comfort.
3. Close the lid tightly to prevent leakage.
4. Insert the nozzle into the irrigator head until it clicks into place.
5. Make sure it is firmly in place.

Nozzle Types

Standard nozzle

Use: General daily cleaning between teeth and gums..

Ideal for: Daily use in people without special needs

Tongue Cleaner Nozzle (Tongue Scraper)

Use: Cleaning the surface of the tongue to remove bacteria that cause bad breath (halitosis).

How to use: With the irrigator in gentle mode, it is passed over the tongue from back to front. Do not press excessively.

Periodontal Tip

Use: Gentle cleaning on periodontal pockets and under the gums..

Ideal for: People with sensitive gums or periodontal disease. It usually has a thinner or softer tip so as not to damage the tissue.

Orthodontic Tip

Use : Cleaning around orthodontic appliances.

Ideal for: Users with brackets or fixed retainers. Tip with specific cleaning channel that removes residues between the devices.

Turn on the device by pressing the power button.

1. Select the desired mode:
 - Soft:** Recommended for sensitive gums or first uses.
 - Normal:** Standard daily use.
 - Pulse:** Deep cleansing with gum massage.

Some models include indicator lights for each mode.

1. Lean slightly over the sink.

2. Place the nozzle in the mouth before turning on the water jet.
3. Close your lips slightly to prevent splashing, allowing water to escape.
4. Direct the jet along the gum line and between the teeth, slowly moving the nozzle.
5. The session usually lasts 45 to 90 seconds.

Tip: Start at the back and work your way up to the front teeth.

When the rinse cycle is complete.

1. Turn off the device after use.
2. Empty the remaining water tank.
3. Remove the nozzle and rinse it with water.
4. Dry the appliance with a soft cloth and store it in a dry place.

WEEKLY CLEANING AND STORAGE

Make sure the device is completely turned off and disconnected from the charger before cleaning it

1. Fill the reservoir with a mixture of warm water and white vinegar (in a 2:1 ratio) or with a mild cleaning solution recommended by the manufacturer.
2. Activate the appliance and let the solution circulate for a few seconds.
3. Empty and rinse the reservoir several times with clean water to remove any residue.
4. Rinse the nozzles with clean water as well. For a deeper clean, you can soak the nozzles in a solution of water and vinegar for 5–10 minutes and then rinse them thoroughly.
5. Replace nozzles every 3 to 6 months or sooner if you notice

wear, deformation, or loss of pressure.

6. Do not use solvents, bleach, alcohol, harsh detergents or hot water above 40°C for cleaning.
7. Do not immerse the irrigator body in water or expose it directly to jets.
8. Do not store the appliance with water in the reservoir and store the irrigator in a dry place, away from sources of heat, direct sunlight and excessive moisture.
9. If you will not be using the appliance for an extended period of time, charge the battery at least once every 3 months to maintain its capacity.

ENVIRONMENTAL PROTECTION – DIRECTIVE 2012/19/EU 2012/19 / UE

In order to preserve our environment and protect health, waste electrical and electronic equipment must be disposed of in accordance with specific standards with the involvement of suppliers and users.

The symbol on the product or packaging indicates that this product should not be disposed of with ordinary household waste.

The user must take the product to a collection point for the recycling of electrical and electronic equipment.

Removal must be in accordance with current local regulations. Consult regional or municipal authorities about options for waste disposal.

TECHNICAL SPECIFICATIONS

Model No.: 5428756

3.7V  1200 mAh / Input 5V  1A

IPX6

This device complies with Directive 2014/35/EU on Low Voltage and Directive 2014/30/EU on Electromagnetic Compatibility.

WARRANTY CONDITIONS

The legal warranty in accordance with current legislation applies to this product. During the legal period, all material and manufacturing defects are covered by the warranty. To make use of the warranty, you must go with the purchase receipt, to the same establishment where you have obtained the item, and choose to refund the money or exchange the defective item for a new one. For after-sales support, you can contact us for any queries in the contact section of <https://adiberia.com/es/contact/questions> Claims for broken accessories or defects due to improper use or care, excessive wear due to intensive non-domestic use, incorrect use due to non-compliance with the instruction manual or normal wear and tear will not be accepted. If this appliance is used incorrectly, for professional or semi-professional purposes or in a manner that is not in accordance with the instructions for use, the warranty is void. Warranty claims will not be accepted if the appliance has been modified or if attempts have been made to repair it. In such a case, any liability in the event of accidents is rejected. Parts subject to wear, as well as cleaning and maintenance or replacement of worn parts are not covered by the warranty, and their replacement is at the buyer's expense. This instruction manual and its future updates can be downloaded digitally at www.adiberia.com

AFTER THE LEGAL WARRANTY PERIOD

Repairs after the warranty period will be at the owner's expense.

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