## SB0166

FR - MASSEUR BOULE MINI - RECHARGEABLE USB

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# Instructions for use

## SB0166 - MINI BALL MASSAGER - USB RECHARGEABLE

### **IMPORTANT, KEEP FOR FUTURE REFERENCE. READ CAREFULLY.**

To ensure safe and proper installation and use, as well as proper functioning of the product, always follow these instructions carefully and keep them for future reference.

### **Safety instructions**

This appliance can be used by children aged 8 years and above and by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they are supervised or have been given instruction concerning the safe use of the appliance and understand the hazards involved. Children should not play with the appliance. Cleaning and user maintenance should not be carried out by children without supervision.

Use the product only as described in this manual. Do not use accessories other than those supplied with the appliance. Before connecting the appliance, ensure that the voltage of the appliance corresponds to the mains voltage. Always switch off the power before any installation work. Do not pull on the power cord to disconnect the appliance.

The appliance has a hot surface. People who are insensitive to heat should be careful when using the appliance.

If the power cord is damaged, it must be replaced by the manufacturer, its after-sales service or similarly qualified persons in order to avoid danger. Never attempt to repair the appliance yourself. Contact a qualified repairer or the shop where you purchased the appliance to carry out any repairs and maintenance.

To prevent the risk of electric shock, do not immerse the device in water or any other liquid.

Do not use the device near water.

Do not insert any objects into the device.

Unplug the appliance and wait for it to cool completely before servicing or cleaning it.

For indoor use only.

This appliance contains batteries that are not replaceable.

The device must only be powered by the very low safety voltage corresponding to the marking on the device.

### **Package contents:**

- 1 massage device • 1 charging cable • 1 user manual

### **Description:**

#### **1. ON/OFF - MULTIMODE**

Speeds: 1-Slow, 2-Medium, 3-Fast, 4-Intermittent.

#### **2. Power input/charging cable (USB-C)/battery (not included)**

### **Features:**

Dimensions: 11.5 x 5 x 10 cm / Battery: 3.7V - 500 mAh / Charging time: 2 hours / Battery life: 1 hour 40 minutes / USB-C charging cable included

### **Charging instructions:**

Connect the Type C charging cable to the device's charging port.

Plug the other end of the charging cable into a rechargeable USB port or connect the AC adapter, then plug it into a power outlet to charge the device.

Charging time is approximately 2 hours. Please keep the charging port dry.

Do not use the device while charging. After charging, unplug the power adapter and start using the device.

Press the power button to start.

**Use:**

This massage gun relieves tension in the body through vibrations and pulses. It is ideal after an intense workout.

Note: to protect your joints and muscles, avoid massaging the same area for too long.

- Press button "1" to turn on the device. The "slow" mode is activated.
- Press the "1" button a second time to switch to "medium" mode.
- Press button "1" a third time to switch to "fast" mode
- Press the "1" button a fourth time to switch to "intermittent" mode
- Press one last time to switch off the device.

**Cleaning and maintenance:**

The device can be cleaned with a clean, dry cloth.

Do not immerse the product in water.

Do not use detergent as this may damage the product.

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