EN User Manual

Cavitation Peeling Device
Model MPK

Thank you for purchasing our product. Before first use, please carefully read this user manual and follow the instructions contained herein. Keep this manual for future reference, as it contains important information.

Package Contents

1. Cavitation device
2. USB-A to USB-C power cable
3. User manual

Check whether the product is damaged. Damaged equipment may be returned to the point of purchase for repair or replacement in case of defects or damages caused by factors inherent in the product itself. In consideration of environmental protection, please dispose of the product packaging in appropriate selective waste containers. Packaging should be kept out of the reach of children. For children's safety, do not leave packaging parts (plastic bags, cartons, styrofoam, etc.) freely accessible. Risk of suffocation!

Intended Use of the Product

The device is intended for skin cleansing, accelerating its nourishment and regeneration. Cavitation is a physical phenomenon involving the formation of microscopic gas bubbles and their sudden collapse, accompanied by pressure changes that result in breaking down dead cells of the horny layer of the epidermis. This occurs when ultrasonic waves are applied to moist skin. The product is intended for home use indoors only and must not be used for professional purposes.

Safety Instructions

- Before use, carefully read the instructions provided. Failure to comply with these instructions may be dangerous or illegal.
- Use the device only as intended.
- The product is designed for facial and neck skin care. It is not intended for laboratory or industrial use. It is not a medical device.
- Packaging is not a toy. Keep it away from children immediately after unpacking.
- The device is not a toy; ensure children do not play with it. The product is intended for adult use only.
- Repairs may only be performed by a qualified manufacturer's service or authorized service center. Repairs by unauthorized or unqualified services may damage the device and void the warranty.
- The manufacturer is not responsible for damages caused by unauthorized modifications.
- Do not use the product if it is damaged.
- Protect the device from drops and strong shocks.
- Prevent water from entering the device. Do not immerse it in water or other liquids.
- Do not use near sinks, bathtubs, or other containers filled with water.
- If any allergic skin reaction occurs during or after use, discontinue use and consult a doctor.
- Use the device at room temperature conditions.
- Wash and dry your hands before use.
- Do not use multiple similar devices simultaneously.
- Warning! Do not use the product more than 3 times per week, and no longer than 3 minutes per session.
- The device uses ultrasound, i.e., high-frequency vibrations. This may produce a mild resonance effect on the skin with low-energy, high-frequency vibrations, enhancing cell vitality, deepening sweat glands, and reducing pores.
- The device's effect may vary depending on skin type, usage range, and frequency.
- Before use, consultation with a specialist is recommended to check for contraindications.
- Warning! Do not press the device too hard against the skin, do not hold the spatula (4) in one place on the skin for too long; move the spatula slowly.
- While charging, connect the device so that the cable cannot be pinched, does not hang loosely, and does not rest on sharp edges.
- Connect the device to an easily accessible outlet to allow quick disconnection in case of emergency.
- The power adapter is not included. Connect the power cable only to sources matching the specifications.
- If unusual noises or smoke appear, immediately disconnect the device from power.
- Clean the device only with a dry, soft cloth. Do not use highly acidic or alkaline cleaning agents.

Contraindications


- Do not use on irritated skin, wounds, active acne, or infections. Do not use with strong suntan.
- Fever conditions.
- Do not use if you have congenital heart defects or circulatory failure.
- Severe general condition, blood circulation disorders.
- Do not use during pregnancy.
- Cancer or post-operative states after tumor removal.
- Presence of metallic foreign bodies in tissues, e.g., dental braces or pacemakers.
- Do not use if you have facial nerve disorders.
- Do not use in case of viral, bacterial, or fungal skin diseases.
- Do not use immediately after tanning beds, sun exposure, or sauna.
- Do not use after alcohol or drug consumption, or if excessively tired or physically weakened.

Battery Safety Instructions

- The device is powered by a built-in rechargeable battery, which is integral and not user-replaceable.
- Avoid exposing the battery to very low or very high temperatures (below 10°C / 50°F or above 45°C / 113°F). Extreme temperatures may affect battery life.
- Avoid contact of the battery with liquids or metal objects, as this may cause total or partial damage.
- Use the battery only for its intended purpose. Do not destroy, damage, or throw into fire — this is dangerous and may cause fire or explosion hazard!
- Swallowing the battery can be fatal! Keep away from children and pets. Seek immediate medical help if swallowed.
- Do not dispose of used or damaged batteries with household waste. Deliver them to designated selective collection points in your community.
- Do not dismantle, modify, open, or short-circuit battery terminals!

- *Note*: If not using the device for a long time, charge it at least once every 3 months. Avoid complete discharge.
- *Note*: The device battery is non-replaceable and an integral part of the device. Do not attempt to replace it.
- Before first use, charge the battery according to the user manual instructions.

Device Components

1. Power button
2. Display
3. USB-C power socket
4. Spatula
5. Spatula cover

Specifications

Power supply:	5 V⎓ 1 A
Battery:	Li-ion 3.7 V⎓ 400 mAh 1.48 Wh
Charging port:	USB-C
Charging time:	approx. 1.5 hours
Operating time:	approx. 60 minutes
Dimensions:	16,7 × 4,1 × 1,7 cm
Weight:	64 g

Charging the Device Battery

Before first use, charge the device battery. Connect the USB-A plug of the power cable to a power source matching the specifications (5 V⎓). Insert the USB-C plug into the device power socket (3). Then plug the power adapter into the mains socket. The battery icon on the display (2) will pulse. When fully charged, the icon will stop pulsing. Disconnect the device from power. Low battery is indicated by the battery icon pulsing in red.

NOTE! The charger is not included.

Usage

To turn the device on or off, press and hold the power button (1). To change the operating mode, press the power button (1) repeatedly. The active function is indicated on the display (2). The device will automatically switch off after 10 minutes of inactivity.


NOTE: It is recommended to test the device on the wrist before using it on the face.

Device Operation


1. Wash your face and ensure it is well-moistened.
2. Turn on the device using the power button (1).
3. Press the button (1) repeatedly to change the mode.
4. Perform cleansing according to the illustrations below, corresponding to the device's functions. Start from the nose, moving from bottom to top. Then cleanse the area from the jawline to below the eyes and cheeks, from bottom to top on both sides. Next, the chin, area above the lips, and forehead, moving the spatula from bottom to top. Do not use directly near the eyes.

Device Functions


1. **Deep Cleansing Mode (Peel):** The spatula produces vibrations to gently remove blackheads, keratin, acne, dead skin cells, and makeup residues. When using Peel mode, a cleansing gel can be applied to create foam to be removed with the device.


2. **Skin Firming Mode (Lift):** Supports blood circulation and accelerates skin metabolism to make it firm, smooth, and delicate. This function can be used with firming serum or cream.


3. **Positive Ion Cleansing Mode (Ion+):** Removes stubborn blackheads, cleanses pores, exfoliates dead and dull skin cells. Use with toner or cleansing milk.

4. **Negative Ion Infusion Mode (Ion-):** Helps achieve elastic skin and reduce wrinkles. Works more effectively with skincare products such as cleansing foam, toner, cream, or cleansing oil. After thorough makeup removal, apply serum, vitamin ampoules, or cream intended for such treatments. Apply the device to prepared skin and start cleansing.

Treatment According to Face Map

Maximum effect with no risk of errors.

The illustration shows precise spatula movement directions during the treatment. Following markings (1–8) helps deeply cleanse pores, stimulate microcirculation, enhance facial contour, and smooth fine wrinkles. Move the spatula according to the arrows at about a 45° angle. This positioning strengthens the device's effect and provides an even lifting effect. It's a great addition to daily anti-aging care—you can act like a professional without leaving home.

Maintenance and Storage

To clean the device, gently wipe it with a dry, soft cloth. Do not use corrosive chemicals or immerse the device in water or any other solution. Follow general safety rules and exercise caution during cleaning. Store the device in a dry, ventilated place away from heat sources. If you do not intend to use the device for a prolonged period, keep it in the original packaging.

Disposal

The device is marked with a crossed-out wheeled bin symbol in accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). Products marked with this symbol must not be disposed of as unsorted municipal waste after the end of their service life. Users are obliged to dispose of used electrical and electronic equipment at designated collection points where such hazardous waste undergoes recycling processes. Collecting such waste in designated places and proper recovery procedures help protect natural resources. Proper recycling of electrical and electronic equipment benefits human health and the environment. For information on where and how to safely dispose of used electronic equipment, contact local authorities, waste collection points, or the place of purchase.

According to Regulation (EU) 2023/1542 on batteries and waste batteries, this product is marked with a crossed-out waste bin symbol indicating that batteries used in this product must not be disposed of with normal household waste but handled according to directives and local regulations. Batteries must not be discarded with unsorted municipal waste. Battery users must utilize the available collection network for battery return, recycling, and disposal. Battery collection and recycling in the EU are subject to separate procedures. For information on battery recycling procedures in your area, contact local municipal offices, waste management institutions, or landfills.

Symbol Explanations

The product complies with European technical standards known as Harmonized European Standards. It complies with Directive 2014/30/EU. A declaration of conformity has been issued, and the CE mark is affixed.

The symbol indicates that the product must not be disposed of with other household waste. Electrical equipment that is not recycled pollutes air, soil, and groundwater. Polish law prohibits mixing electrical and electronic waste with other waste under penalty of fine.

The symbol indicates the type of material used for packaging — flat cardboard.

The symbol indicates that the packaging is recyclable.

The symbol indicates that the user should read the instruction manual before use.

Manufacturer

TelForceOne S.A.
Krakowska 119, 50-428 Wrocław, Poland
Hotline: (+48 71) 71 77 400, e-mail: info@telforceone.com, www.telforceone.pl
Made in China

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